ID **Italian**
design

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EN PRECAUTIONS

Please read the instructions carefully before using the appliance.

To reduce any risk of accident, take the following precautions:

1. Never use the appliance while you are in the bath.

2. IEC 60417

- 5582: Do not

place or store the

appliance anywhere where it could fall or be pushed towards the bath or sink.

3. Do not immerse or drop into the water or any other liquid.

4. If the appliance falls into the water, unplug it immediately. Do not place your hand in the water.

5. Never use with wet hands.

6. Keep the appliance away from water.

7. An appliance should never be left unattended while it is plugged into the mains.

8. This appliance can be used by children aged from 8 years and

above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.

9. Use this appliance only for its intended purpose, as indicated in this manual. Do not use any accessory that is not recommended by the manufacturer.

10. Never use this appliance if the power cord or plug is damaged, if it is malfunctioning, if it has been dropped or has fallen into water. Take the appliance to your nearest service centre for inspection or repair.

11. Keep the power cord away from sources of heat or hot surfaces. Do not coil the power cord around the appliance.

12. Never drop or insert any object into the openings.

- 13.** Do not manipulate the appliance.
- 14.** This product must not be used by children or by other people without assistance or supervision if their physical, sensory or mental condition limits safe use of the appliance.
- 15.** This appliance does not contain any user-serviceable parts. Do not try to repair the appliance yourself. Always take it to a qualified technician.
- 16.** Do not use this appliance while walking or driving.
- 17.** Do not use this appliance in a hot and humid place.
- 18.** Do not use this appliance if you have symptoms of eye disease of any kind or eye surgery. Its use is strictly prohibited if you suffer from glaucoma, cataracts, or retinal detachment.
- 19.** This device should be used with caution if you have had a stroke, hypertension, or hypotension. Consult your doctor before use.
- 20.** If you feel any discomfort or pain during use, please stop using it immediately.
- 21.** People with physical or mental discomfort or mental handicap should use this device with caution. Use it under the supervision of accompanying personnel. If you feel any discomfort during use, please stop using it immediately.
- 22.** Avoid using the device if you have been drinking alcohol or experiencing physical discomfort.
- 23.** This appliance heats up slightly during use (42°C) and should be used with care when used with persons sensitive to temperature.
- 24.** Put the appliance on before turning it on to avoid damaging the airbags.
- 25. IMPORTANT:** Before using make sure you have removed contact lenses, hairpins, hair accessories, etc.

26. This massage device is a relaxation product and cannot replace drugs or medical equipment.

27. This device is designed as an electronic care product for personal use as a massager and not for use in diagnostic, therapeutic, and other commercial applications.

28. Clean this appliance with a cotton cloth soaked in medical alcohol. Do not wash the product with water.

29. Use only the supplied accessory for your charge.

30. Use this appliance with your eyes closed. After use, your vision may be blurred. After 1-3 minutes your eyesight will be normal again.

EN BEFORE BEING USED

Thanks to its high technology, the ocular pressotherapy massager combines heat, massage, and acupressure techniques on the eye and temples, emitting gentle vibrations that improve blood circulation and that together with the sound effects emitted by its speakers achieve an

immersion effect for a total physical and mental relaxation.

Acupressure is done by means of air pressure so it is a very gentle and delicate pressure with the eyes. This technique allows you to relax the eye and head nervous system since much of the tension is concentrated in the temples.

The massager works at a pleasant constant temperature of 42° which, in addition to helping relaxation, also helps to activate blood circulation, improve the ability to supply nutrients and help relieve dry eyes.

The different sound effects act for 10 minutes and facilitate sensory immersion and enhance the relaxing effects of the massage.

This device is specially designed for people who, in one way or another, strain their eyes for a long time to be in front of a screen, text, driving, wear glasses, etc, whether for studies, work or leisure. They are also recommended for insomniacs.

- Make sure the appliance is charged before use, otherwise it will not work. (with low power). Lo puesto entre parentesis, sobra.
- Do not use this appliance if the low battery alarm (voice) sounds.
- Make sure the appliance is switched off during the charging process. This appliance has a self-protection program by which it will stop automatically if it is in use.
- This device has a built-in lithium battery whose charging time is 2 - 3 hours.

- Please connect the adapter with the USB port.
- The red light indicates that the device is charging.
- The blue light indicates that the device is already charged.
- Be sure to unplug the adapter after charging to avoid charging for a long time and shortening the battery life.
- This device can be charged by the power supplied by a computer or a telephone or an external battery, although it is recommended that it be carried out using a standard adapter plugged into the mains.

EN HOW TO USE

1. Each session lasts 10 minutes, after which the device will turn off automatically. If you wish to extend the massage, you must (turn the appliance back on.) switch on the appliance again.

1. To turn on the appliance, press and hold the on/off button for 3 seconds. The device will start by default in automatic mode. Pressing the on / off button again (will select) the next mode (dynamic mode) will be selected and pressing again (will select) the last mode (rest mode) will be the one to work.

A. Automatic mode: Air pressure + heating + river sound

B. Dynamic mode: Air pressure + sound waves

C. Rest mode: Heating + sound birds

EN ENVIRONMENT

Disposal of old Electrical & Electronic Equipment (Applicable throughout the European Union and other European countries

with separate collection programs. This symbol, found on your product or on its packaging, indicates that this product should not be treated as household waste when you wish to dispose of it. Instead, it should be handed over to an applicable collection point for the recycling of electrical and electronic equipment. By ensuring this product is disposed of correctly, you will help prevent negative consequences to the environment and human health, which could otherwise be caused by inappropriate disposal of this product. The recycling of materials will help to conserve natural resources. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

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